POLYVIEW® PLUS HRP (anti-mouse) reagent

High sensitivity, low background nanopolymer detection reagent for use with HIGHDEF® chromogens in ISH and IHC applications

The POLYVIEW® PLUS HRP (anti-mouse) reagent is a ready-to-use, biotin-free, one-step detection reagent suitable for immunohistochemical detection of antigens in formalin-fixed paraffin-embedded tissues and frozen sections. The detection reagent may also be used with blood smears, cytosmears, and cell preparations. In addition, the POLYVIEW® PLUS HRP (anti-mouse) reagent can be adapted for *in situ* hybridization detection.


Citations: 6 View Online »

Ordering Information Order Online »


ENZ-ACC104-0150	150 tests

Manuals, SDS & CofA View Online »


- Ready-to-use reagent
- Biotin-free nanopolymer detection circumvents endogenous biotin background
- High intensity color development delivers sharp, crisp staining
- Compatible with both automated and manual platforms
- The Antibody Blocker/Diluent further reduces background (sold separately)
- Customize with HIGHDEF® chromogens of your choice (sold separately)




Dual IHC staining performed on formalin-fixed paraffin-embedded tissue section from human brain with Alzheimer's disease. The tissue sections were stained for alpha-amyloid with POLYVIEW® PLUS AP (anti-rabbit) reagent (Prod. No. ENZ-ACC110) and HIGHDEF® Green AP Chromogen/Substrate (Prod. No. ENZ-ACC130), and phospho-tau with POLYVIEW® PLUS HRP (anti-mouse) reagent (Prod. No. ENZ-ACC104) and HIGHDEF DAB Chromogen/Substrate Set (Prod. No. ENZ-ACC105). Picture courtesy of Dr. Karelle Leroy Laboratoire d'histologie, de neuroanatomie et de neuropathologie, Université Libre de Bruxelles, Bruxelles, Belgium.




Tonsil tissue stained with anti-Ki-67 antibody (mouse), followed by POLYVIEW® PLUS HRP Reagent






Cervical tissue stained with anti-Ki-67 antibody (mouse), followed by POLYVIEW® PLUS HRP (anti-mouse) reagent.



Dual IHC staining performed on formalin-fixed paraffin-embedded tissue section from transgenic mouse brain expressing human mutated APP and PS1. The tissue sections were stained for alpha-amyloid with POLYVIEW® PLUS AP (anti-rabbit) reagent (Prod. No. ENZ-ACC110) and HIGHDEF® Green AP Chromogen/Substrate (Prod. No. ENZ-ACC130), and phospho-tau with POLYVIEW® PLUS HRP (anti-mouse) reagent (Prod. No. ENZ-ACC104) and HIGHDEF DAB Chromogen/Substrate Set (Prod. No. ENZ-ACC105). Picture courtesy of Dr. Karelle Leroy Laboratoire d'histologie, de neuroanatomie et de neuropathologie, Université Libre de Bruxelles, Bruxelles, Belgium.



Cervical tissue stained with anti-Ki-67 antibody (mouse), followed by POLYVIEW® PLUS HRP (anti-mouse) reagent

Handling & Storage

Use/Stability	As indicated on product label or CoA when stored as recommended.
Handling	Do not freeze.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application

IHC, ISH (in situ hybridization)

Technical Info / Product Notes

The POLYVIEW® PLUS HRP (anti-mouse) reagent has been developed by directly labeling anti-mouse immunoglobulins with a nanopolymer of horseradish peroxidase (HRP) using proprietary technology. This ensures consistent and reproducible immunodetection of mouse primary antibodies against nuclear, cytoplasmic and membrane antigens in different types of tissues. The single step methodology enables faster staining procedures than traditional two-step methods using biotin and avidin/streptavidin conjugates, with significantly lower background.

The POLYVIEW® PLUS HRP (anti-mouse) reagent is suitable for use with mouse IgG antibodies. The reagent can be used for manual staining or with automated staining instruments and are well suited for multiplex immunohistochemical staining assays.

The Antibody Blocker/Diluent (ENZ-ACC108) further reduces background and is suitable for diluting the POLYVIEW® PLUS HRP (anti-mouse) reagent.

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.