PROTEOSTAT®
Protein aggregation assay

Quantitative Detection of Protein Aggregates from Visible to Subvisible Particles

PROTEOSTAT® Protein aggregation assay provides a simple, homogenous assay format for monitoring peptide and protein aggregation in solution. This is useful for defining optimal storage formulations for proteins, for screening of compounds that promote or inhibit protein aggregation and potentially for the sensitive measurement of molecular chaperone activity. The assay can be employed to streamline protein processing and formulation optimization procedures. Relative to conventional protein aggregation detection dyes, such as Thioflavin T, the Enzo PROTEOSTAT® detection reagent can detect aggregates from a broader range of proteins, yields a much brighter signal, provides at least 2 orders of magnitude linear dynamic range and offers superior performance across a broad range of pH values (4~10) and buffer compositions. Sensitivity for this assay is in the sub-micromolar range and as little as 1-5% protein aggregate is detectable in a concentrated protein solution. The assay is capable of providing quantitative analysis of protein aggregation in a robust and high-throughput fashion (Z' factor score >0.5). Lyophilized native and aggregated protein are provided as negative and positive controls for monitoring changes in protein aggregation status.

PROTEOSTAT® Protein Aggregation Standards kit (Prod. No. ENZ-51039) is the *only commercially available* protein aggregation standards assay with stabilized, high-quality reference samples for generating trace protein aggregate levels in concentrated monomeric IgG.
Easy to use - *simply add water!*

Citations: 85 View Online »

- A simple, sensitive, homogenous fluorescent assay
- Validated for use with microplate or flow cytometry platform
- Extensively benchmarked with IgG
- Optimize buffers and excipients for protein formulation
- Performs with a wide pH and ionic strength range
- Use with PROTEOSTAT® Protein Aggregation Standards for accurate quantification of aggregated protein in solution.

Ordering Information

Order Online »

ENZ-51023-KP050	50 tests
ENZ-51023-KP002	2x96 tests

Manuals, SDS & CofA

View Online »







Figure 2. Flow Cytometry Application: PROTEOSTAT® Dye Mixed with Bodipy Dye (Pyrromethene 546) Differentiates Oil Droplets From True Protein Aggregates. Bodipy vs PROTEOSTAT® fluorescence of IgG aggregates (A), Silicon oil droplets (B) and a mixture of Silicon oil droplets and IgG aggregates (C).




Effective linear dynamic range for antibody aggregate detection using PROTEOSTAT® Detection Reagent compared with Thioflavin T. Relative fluorescence unit values (RFUs) may differ depending upon the microplate reader employed for the analysis.





Dye is immobilized when bound to the aggregate and begins to fluoresce.



Thioflavin T: early prototype dye in the design of PROTEOSTAT® assay which also rotates around a single bond (red arrow) in the absence of protein aggregates.


Lysozyme
Lysozyme + tri-N-acetylchitotriose
RFU
Hours
25,000
20,000
15,000
10,000
5,000
0
0
100
200
300
400
500
600
700
800


Aggregate Signal
Time
Lag Phase
Growth Phase
Plateau Phase
Region is considered
100% aggregated.
0 hrs
6 hrs
12 hrs
18 hrs
24 hrs


1
0.9
0.8
0.7
0.6
0.5
0.4
0.3
0.2
0.1
0
400
450
500
550
600
650
700

Handling & Storage

Use/Stability	With proper storage, the kit components are stable up to the date noted on the product label. Store kit at -20°C in a non-frost free freezer or -80°C for longer term storage.
Short Term Storage	-20°C
Long Term Storage	-80°C
Shipping	Dry Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	Flow Cytometry, Fluorescence microscopy, Fluorescent detection, Microplate
Application Notes	This kit has been specifically designed for monitoring of protein aggregate formation in solution.
Contents	PROTEOSTAT® detection reagent PROTEOSTAT® positive control PROTEOSTAT® negative control 10X PROTEOSTAT® assay buffer
Quality Control	A sample of PROTEOSTAT® Protein aggregation assay was used to assay (1) 20 µM aggregated lysozyme, (2) 20 µM native lysozyme; and (3) mixtures of 5% aggregate. The Z' factor is greater than 0.5 and the 5% aggregate signal is greater than 3 standard deviations above the no-aggregate control.
Quantity	For ENZ-51023-KP050: 50 tests in a 96 well plate or 16 flow cytometry tests For ENZ-51023-KP002: 2 x 96-well tests or 70 flow cytometry tests

Technical Info / Product Notes

Application Notes:

Particle analysis of therapeutic protein formulations with ImageStreamX® Imaging Flow Cytometry and the PROTEOSTAT® Protein Aggregation Assay

Prediction of Aggregation Propensity and Monitoring of Aggregation of Antibody-Drug Conjugates (ADC) using ProteoStat® Reagents

Enzo and PROTEOSTAT are trademarks of Enzo Life Sciences, Inc. Several of Enzo's products and product applications are covered by US and foreign patents and patents pending.

Cited samples:

PROTEOSTAT® Cited Samples

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.