2′-AHC-cAMP-Agarose

Cyclic AMP immobilized on agarose gel

The protein kinase A activator cAMP immobilized on agarose by an aminohexylcarbamoyl spacer attached to the ribose 2' position of the ligand. The gel can be used for affinity chromatography of various cyclic nucleotide-responsive proteins, especially those which tolerate modification of the 2' ribose hydroxyl group, such as some phosphodiesterases. This type of gel is also available with immobilized PDE-resistant Rp-cAMPS (Rp-2'-AHC-cAMPS-Agarose (Prod. No. BLG-A049)) or Sp-cAMPS (Sp-2'-AHC-cAMPS-Agarose (Prod. No. BLG-A050)).


Ordering Information

Order Online »

BLG-A058-25	2.5ml
BLG-A058-60	6ml

Manuals, SDS & CofA

View Online »


NH2
N
N
N
N
O
O
O=P—O
O
O
H
N
N
H
O
Na+
O−

Handling & Storage

Long Term Storage	+4°C
Shipping	Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	2′-(6-Aminohexylcarbamoyl)adenosine-3′,5′-cyclic monophosphate; immobilized on agarose gel
Couple Type	Modified nucleotides
Technical Info / Product Notes	For the Original Manufacturer’s data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.