2′-AHC-cUMP-Agarose

Cyclic UMP immobilized on agarose gel

The putative second messenger cyclic UMP immobilized on agarose by an aminoethylcarbamoyl spacer attached to the ribose 2'-position of the ligand. The gel can be used for affinity chromatography of cUMP-responsive proteins, such as protein kinases, phosphodiesterases and others. Detailed technical information available.


Ordering Information

Order Online »

BLG-A164-06	0.6ml
BLG-A164-25	2.5ml
BLG-A164-60	6ml

Manuals, SDS & CofA

View Online »


O
HN
O
N
O
O
O=P—O
O
Na+
O−
O
H
N
O
N
H

Handling & Storage

Long Term Storage	+4°C
Shipping	Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	2′-O-(6-Aminohexylcarbamoyl)uridine-3′,5′-cyclic monophosphate, immobilized on agarose
Couple Type	Modified nucleotides
Formulation	Suspension in 30 mM Na_2HPO_4 buffer (pH 7).
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.