# 2′-O-Me-c-diAMP

## Analogue of c-diAMP

*BLG-M085-005 (5 x 0.1 µmol pack size) is not sold in the U.S. or Canada. Please contact us for available options.*


## Ordering Information

Order Online »

| | |
|---|---|
| **BLG-M085-001** | 0.1µmol |
| **BLG-M085-005** | 5x0.1µmol |

## Manuals, SDS & CofA

View Online »


NH2
N
N
A
N
N
Na+
O−
O=P—O
5'
O
O
O
3'
3'
O
OH
O
5'
O—P=O
O−
Na+
N
N
A
N
N
NH2

## Handling & Storage

| | |
|---|---|
| **Long Term Storage** | -20°C |
| **Shipping** | Ambient Temperature |

## Regulatory Status

RUO - Research Use Only

## Product Details

| | |
|---|---|
| **Alternative Name** | 2′-O-Methyl-cyclic diadenosine monophosphate . sodium salt |
| **Formula** | $C_{21}H_{26}N_{10}O_{12}P_2$ (free acid) |
| **MW** | 672.4 (free acid) |
| **Purity** | > 95% HPLC |
| **Quantity** | 0.1µmol=~67µg |
| **Technical Info / Product Notes** | For the Original Manufacturer's data sheet please click here. |

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.