AMPIVIEW® Pantr1 (AS) Mouse DNP RNA Probes Set

Mouse-specific p53-regulated lncRNA (Pantr1) probe for oncology research

The AMPIVIEW® Pantr1 (AS) Mouse DNP RNA Probes set is a ready-to-use in situ hybridization (ISH) probe, powered by Enzo's LoopRNA ISH™ technology, for chromogenic detection of Pantr1 (Panther 1 / lincRNA-p21), a mouse-specific p53-regulated long non-coding RNA involved in apoptotic signaling and cell cycle arrest pathways. AMPIVIEW® Pantr1 probes enable visualization of this p53-downstream lncRNA in murine tissue for cancer biology research and preclinical oncology studies. in formalin-fixed paraffin-embedded (FFPE) tissue sections and cells. The probe is DNP-labeled and detected using an unconjugated anti-DNP linker followed by a polymer-based POLYVIEW® PLUS detection system combined with HIGHDEF® chromogens and counterstain, visualized under a standard light microscope. Supplied as 2 x 1.0 mL ready-to-use probe (ENZ-GEN560) with 1 x 2.0 mL AMPIVIEW® Hybridization Buffer (Fast Hyb Solution) (ENZ-ACC187).

AMPIVIEW® RNA probes are designed to detect the RNA or DNA/RNA of Pantr1 (Panther 1 / lincRNA-p21), a mouse-specific p53-regulated long non-coding RNA involved in apoptotic signaling and cell cycle arrest pathways. AMPIVIEW® Pantr1 probes enable visualization of this p53-downstream lncRNA in murine tissue for cancer biology research and preclinical oncology studies. In situ hybridization (ISH) is a molecular technique that uses labeled nucleic acid probes to detect and localize specific RNA or DNA sequences directly within tissue sections or cells. Powered by Enzo's LoopRNA ISH™ technology, AMPIVIEW® RNA probes deliver superior sensitivity for chromogenic detection of RNA, DNA, or RNA/DNA targets on formalin-fixed paraffin-embedded (FFPE) tissue and cells, and are compatible with Enzo's complete immunohistochemistry detection workflow, from antigen retrieval to counterstain.

- Powered by Enzo's LoopRNA ISH™ technology for superior sensitivity
- Ready-to-use, DNP-labeled probe compatible with AMPIVIEW® Rabbit Anti-DNP Linker, POLYVIEW® PLUS detection and HIGHDEF® chromogens
- 1-hour hybridization with AMPIVIEW® Hybridization Buffer (Fast Hyb Solution) (ENZ-ACC187) for fast, same-day results
- Compatible with manual staining and automated platforms (LEICA BOND, Roche DISCOVERY/ULTRA) with prior optimization

Ordering Information

Order Online »

ENZ-GEN560-2000	

Manuals, SDS & CofA

View Online »

Handling & Storage

Use/Stability	Stable under recommended storage conditions up to the expiration date indicated on the vial label. Avoid freeze-thaw cycles.
Handling	Protect from light and humidity. Avoid freeze-thaw cycles. Avoid RNases; wear gloves and work in a clean environment.
Short Term Storage	-20°C
Long Term Storage	-20°C
Shipping	Dry Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	ISH
Contents	2 x 1.0 mL AMPIVIEW® Pantr1 (AS) Mouse DNP RNA Probes set (ENZ-GEN560); 1 x 2.0 mL AMPIVIEW® Hybridization Buffer (ENZ-ACC187)
Quantity	2 x 1.0 mL probe + 1 x 2.0 mL hybridization buffer
Recommendation Dilutions/Conditions	Ready to use solution for manual use; Optimization may be required for automated staining systems.
Specificity	Pantr1 / lincRNA-p21 (antisense)
Technical Info / Product Notes	Ready-to-use probe for manual staining. Prior optimization required for automated staining systems. AMPIVIEW® RNA probes have been optimized with Enzo's AMPIVIEW® and POLYVIEW® Detection Systems combined with HIGHDEF® chromogens. Use of other linkers or polymer-based reagents will require prior optimization. Autostainers such as LEICA BOND or Roche DISCOVERY/ULTRA can be used but may require prior optimization

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.