Antigen Retrieval Reagent, pH 9 (10X)

Antigen Retrieval Reagent, pH 9 (10X) is used to recover antigens masked by fixation in cross linking fixatives such as formalin.

Citations: 1 View Online »

Ordering Information

Order Online »

ENZ-ACC113-0100	100ml

Manuals, SDS & CofA

View Online »



Antigen Retrieval Reagent, pH 9 (10X)

Handling & Storage

Use/Stability	As indicated on product label or CoA when stored as recommended.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Formulation

Liquid. In Tris-EDTA buffer.

Technical Info / Product Notes

In order to perform immunostaining, tissue specimens should be preserved in an appropriate fixative. Fixation stops tissue autolysis, preserves tissue structures, and immobilizes antigens. However, antigens undergo chemical alteration of their primary, secondary and tertiary structures during fixation. Antigenic sites may be masked due to changes induced in the epitope or neighboring proteins. Enzymatic treatment with proteolytic enzymes (i.e. pepsin, trypsin or pronase) has been performed to expose the masked antigens.

Shi et al. (1991) have reported that treating tissue sections with a heavy metal solution in a microwave oven can cover masked antigens significantly. However, heavy metals in the solution increase the risk of toxic exposure to lab personnel. To reduce the risk, we have developed an antigen unmasking solution which is free of heavy metals. The reagent is provided in a convenient 10X concentrate. Use of this antigen recovery buffer avoids unnecessary heavy metal exposure and other handling and disposal issues.

Interpretation of the results will be the sole responsibility of the user.

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.