Biotin16 NT Labeling Kit

Biotin16 Nick Translation Labeling Kit contains all reagents (except template and materials for purification of the probe) required for nick translation-based Biotin labeling of DNA providing a highly efficient, easy-to-perform and rapid labeling technology. The labeling principle is similar to the Biotin-Nick Translation Mix (Roche). Nick translation labeling is based on the reverse activities of Polymerase I and DNase I. DNase I is able to introduce randomly distributed nicks to double stranded DNA at low enzyme concentrations. The 5'→3' exonuclease activitiy of Polymerase I removes nucleotides from the 3' side of the nick while synthesizing a partial new complementary strand using the 3'-OH termini as primer. In the presence of labeled dUTP Polymerase I incorporates labeled dUTP instead of dTTP. The BIO16 NT labeling mix contains an optimized Biotin-labeled dUTP for incorporation into DNA by nick translation using DNA Polymerase I. Biotin-16-dUTP is enzymatically incorporated into DNA as substitute for its natural counterpart dTTP. Optimal substrate properties and thus labeling efficiency as well as an efficient detection of the Biotin moiety is ensured by a 16-atom linker attached to the C5 position of uridine. The well balanced polymerase / nuclease activities of the enzyme mix ensure the generation of highly labeled double stranded DNA fragments. The resultant Biotin-labeled DNA is suitable for application in FISH and other nucleic acid hybridization assays.

Ordering Information

Order Online »

JBS-PP-310S-BIO16	10Reactions

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	-20°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Quantity	20 µl each
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.