C8-Dihydroceramide

Negative control for C8 ceramide


C_8 Dihydroceramide is inactive and may be used as a negative control for C_8 ceramide.

Citations: 1 View Online »

Ordering Information Order Online »

BML-SL113-0025	25mg

Manuals, SDS & CofA View Online »


O
HN
CH3
HO
CH3
OH

Handling & Storage

Use/Stability	As indicated on product label or CoA when stored as recommended.
Long Term Storage	-20°C
Shipping	Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Appearance	White waxy solid.
CAS	145774-33-0
Formula	$C_{26}H_{53}NO_3$
MW	427.7
Purity	≥98% (TLC)
Solubility	Soluble in DMF (25 mg/ml), DMSO (5 mg/ml) or ethanol (10mg/ml).

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.