Direct 8-iso-$PGF_{2\alpha}$ ELISA kit

Highly sensititive ELISA, measuring 8-iso-Prostaglandin $F_{2\alpha}$ in plasma, serum and tissue.

The measurement of F2-isoprostanes has emerged as one of the most reliable approaches to assess oxidative stress status *in vivo*, providing an important tool to explore the role of oxidative stress in the pathogenesis of human disease. The is a colorimetric competitive immunoassay kit for the quantitative determination of 8-iso-Prostaglandin $F_{2\alpha}$ in biological fluids such as plasma, serum and tissue. To measure isoprostane in urine or culture supernatants, use the 8-iso-$PGF_{2\alpha}$ EIA Kit.

8-iso-$PGF_{2\alpha}$ is produced *in vivo* by both non-cyclooxygenase and cyclooxygenase dependant mechanisms from arachidonic acid. It has been shown to be a potent vasoconstrictor, a potential mediator of hepatorenal syndrome and a mutagen in 3T3 cells and vascular smooth muscle cells. 8-iso-$PGF_{2\alpha}$ has been shown to circulate in plasma and is excreted in urine.

Prostaglandin $F_{2\alpha}$ ($PGF_{2\alpha}$) is formed in a variety of cells from PGH_2, which itself is synthesized from arachidonic acid by the enzyme prostaglandin synthetase. $PGF_{2\alpha}$ is often viewed as an antagonist to PGE_2 due to their opposing effects on various tissues. $PGF_{2\alpha}$ is a potent bronchoconstrictor and has been implicated in asthma attacks. $PGF_{2\alpha}$ is also involved in reproductive functions including corpus luteum regulation, uterine contractions, and sperm motility. This has led to its use in terminating pregnancies and inducing labor at term. High levels of $PGF_{2\alpha}$ have also been associated with preeclampsia.

Citations: 34 View Online »

- Non-radioactive measurement of isoprostane in serum and plasma
- Flexible options, either overnight or same day results
- High sensitivity, measuring as little as 40 pg/ml
- High throughput format with results for up to 39 samples in duplicate in less than 3 hours


Ordering Information

Order Online »

ADI-900-091	96 wells

Manuals, SDS & CofA

View Online »



Direct 8-iso-PGF2α ELISA Kit Graph for <3 hour format

Handling & Storage

Long Term Storage +4°C

Shipping Blue Ice

Regulatory Status RUO - Research Use Only

Product Details

Alternative Name	8-iso-Prostaglandin F2α
Application	Colorimetric detection, ELISA
Application Notes	For the quantitative determination of 8-iso-$PGF_{2\alpha}$ in plasma, serum, and tissue from any species. Cited sample types include cell lysate and urine.
Assay Time	<3 hours or Overnight + 45 min
Compatibility	This product is compatible with the **Absorbance 96 Plate Reader**.
Contents	GxR IgG Microtiter plate, Sample diluent, Conjugate, Antibody, Wash buffer concentrate, Standard, pNpp Substrate, Stop solution, Neutralizing reagent
Sensitivity	103.2 pg/mL (<3 hour format), 40.0 pg/mL (overnight) (range 160 – 100,000 pg/mL)
Species Reactivity	Species independent
Technical Info / Product Notes	**Cited samples:** Wildlife Endocrinology: Cited Sample Types
Wavelength	405 nm

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.