IκB kinase substrate (biotinylated)

This peptide sequence, RHDSGLDSMKD, is derived from the area surrounding serine 32 and serine 36 of IκBα, which is phosphorylated by IκB Kinase (IKK). The biotin allows the peptide to be used in IKK kinase assays with streptavidin-bound membranes.

Ordering Information

Order Online »

BML-P148-0001	1mg

Manuals, SDS & CofA

View Online »

Handling & Storage

Use/Stability	As indicated on product label or CoA when stored as recommended.
Long Term Storage	-20°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

MW	1486.9
Sequence	Biotin-Arg-His-Asp-Ser-Gly-Leu-Asp-Ser-Met-Lys-Asp

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.