IHC Wash Buffer (10X)

The IHC wash buffer is a for washing and removing excess reagents during immunostaining.

Optimal immunostaining not only depends on the specificity of the primary antibody and other immunoreagents but also depends on obtaining a good signal to noise ratio. Binding of an antibody to its epitopes involves van der Waals forces, electrostatic forces and hydrophobic forces. Certain antibodies have tendency to bind loosely and nonspecifically to unrelated epitopes, which can create undesired background staining.In order to remove these nonspecifically bound antibodies, a thorough washing is required after each immunostaining step. The IHC wash buffer is specifically designed to remove such loosely bound antibodies effectively and efficiently and to provide a cleaner background staining.

Ordering Information

Order Online »

ADI-950-235-1000	1000ml

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage Ambient

Shipping Ambient Temperature

Regulatory Status RUO - Research Use Only

Product Details

Formulation Liquid. 10X wash buffer.

Technical Info / Product Notes Interpretation of the results will be the sole responsibility of the user.

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.