# Mouse anti-biotin linker (Ready-to-use)

This mouse anti-biotin linker is a ready to use signal amplifying reagent optimized to increase the detection sensitivity of ENZO PATHO-GENE® HPV probes. Although mouse anti-biotin linker was created to enhance ENZO probes, it will recognize any biotin associated nucleotide or macromolecule.

**This antibody is covered by our Worry-Free Guarantee.**

## Ordering Information

Order Online »

| | |
|---|---|
| **ENZ-32893-6000** | 6ml |

## Manuals, SDS & CofA

View Online »

## Handling & Storage

| | |
|---|---|
| **Use/Stability** | As indicated on product label or CoA when stored as recommended. Stable for at least one year after receipt when stored as recommended. |
| **Handling** | Do not freeze. |
| **Long Term Storage** | +4°C |
| **Shipping** | Blue Ice |

## Regulatory Status

RUO - Research Use Only

## Product Details

| | |
|---|---|
| **Application** | ISH (in situ hybridization) |
| **Formulation** | Liquid. Ready-to-use solution in TBS containing 0.05% Kathon and 1% BSA. |
| **Host** | Mouse |
| **Isotype** | IgG1 |
| **Species Reactivity** | Species independent |
| **Worry-free Guarantee** | **This antibody is covered by our Worry-Free Guarantee .** |

**Last modified: May 29, 2024**

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.