Oxytocin ELISA kit

The most sensitive, specific and flexible alternative to more costly and time-consuming methods for oxytocin detection.

The Oxytocin EIA kit is a colorimetric competitive enzyme immunoassay kit with results overnight + 1 hour. Absorbance is read at 405 nm. Very low reactivity with vasopressin provides confidence in assay results. Sensitive quantification of 15pg/mL oxytocin. Ready-to-use liquid color-coded reagents reduce errors in the lab.

Oxytocin is a neurohypophysial peptide which is produced in the paraventricular nuclei of the hypothalamus and stored in the posterior pituitary. The molecule consists of nine amino acids linked with a [1-6] disulfide bond and a semi-flexible carboxyamidated tail. A hormone once thought to be limited to female smooth muscle reproductive physiology, more current findings have determined that oxytocin also functions as a neurotransmitter, may be involved in neuropsychiatric disorders, social/sexual behavior, and is important in male reproductive physiology. Oxytocin and the related neurohypophysial peptide, Arg-Vasopressin, maintain renal water and sodium balance.


Citations: 186 View Online »

Ordering Information Order Online »

ADI-901-153A-0001	5x96 wells

Manuals, SDS & CofA View Online »

- Highly specific ELISA which does not detect vasopressin, providing confidence in assay results
- Ultra-sensitive measurement, detecting as little as 15 pg/mL oxytocin
- Cost effective solution, faster and less costly than LC/MS
- Easy-to-use liquid color-coded reagents and pre-coated plates minimize errors
- Widely published in peer reviewed literature
- Save with bulk kit package!


100
80
60
40
20
0
B/Bo (%)
0.700
0.560
0.420
0.280
0.140
0.000
Optical Density
B/Bo (%)
Optical Density
10
100
1,000
Oxytocin Conc. (pg/mL)

Handling & Storage

Use/Stability	Store all components at +4°C, except standard and conjugate at -20°C.
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	Colorimetric detection, ELISA
Application Notes	For the quantitative determination of Oxytocin in culture supernatants, milk, plasma, and serum from any species. Validated sample types include cerebral spinal fluid, saliva, tissue and urine.
Assay Time	Overnight + 1 hour
Compatibility	This product is compatible with the **Absorbance 96 Plate Reader**.
Contents	GxR Microtiter plate, Conjugate, Antibody, Assay buffer concentrate, Wash buffer concentrate, Standard, pNpp Substrate, Stop solution
Crossreactivity	Mesotocin (7%), Oxytocin (100%), Arg8-Vasotocin (7.5%), and <0.02% for other related molecules.
Sensitivity	15pg/ml (range 15.6 – 1,000 pg/ml)
Species Reactivity	Species independent
Technical Info / Product Notes	**Cited samples:** Hormone ELISAs: Cited Samples Wildlife Endocrinology: Cited Sample Types The Oxytocin ELISA kit has been extensively benchmarked for behavioral studies with plasma samples from a variety of species. **Featured in:** Frontiers in Neuroscience – Validation of a commercially available enzyme immunoassay for the determination of oxytocin in plasma samples from seven domestic animal species.
UniProt ID	P01178 (human)

Wavelength	405 nm

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.