Phi29 DNA Polymerase (L pack)

Bacillus subtilis phage phi29, recombinant, E. coli

Ordering Information

Order Online »

JBS-PCR-381L	5x200U

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	-20°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Appearance	liquid
Application Notes	Rolling-circle amplification (RCA) Multiple displacement amplification (MDA) Whole genome amplification (WGA) Preparation of DNA template for sequencing
Concentration	10 units/µl
Source	Bacillus subtilis phage phi29, recombinant, E. coli
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.