Rp-8-AEA-cAMPS-Agarose

Immobilized cAMP Antagonist

The metabolically stable protein kinase A inhibitor Rp-cAMPS immobilized on agarose by an aminoethylamino spacer attached to position 8 of the ligand. Free Rp-cAMPS binds to the protein kinase A holoenzyme but prevents dissociation to free catalytic and regulatory subunits. In addition, it is not hydrolized by mammalian phosphodiesterases.The gel can be used for affinity chromatography of various cyclic nucleotide-responsive proteins such as protein kinase A (presumably as its holoenzyme), phosphodiesterases and others. This ligand is also available with a longer spacer arm (Rp-8-AHA-cAMPS-Agarose Prod. No. BLG-A012)).


Ordering Information

Order Online »

BLG-A007-06	0.6ml
BLG-A007-25	2.5ml
BLG-A007-60	6ml

Manuals, SDS & CofA

View Online »


NH2
N
N
N
N
H
N
H
N
O
O
O=P
O
OH
S
Na+
−

Handling & Storage

Long Term Storage	+4°C
Shipping	Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	8-(2-Aminoethyl)aminoadenosine-3′,5′-cyclic monophosphorothioate, Rp-isomer; immobilized on agarose gel
Couple Type	Modified nucleotides
Formulation	Liquid. 0.05M Na_2HPO_4 buffer (pH7), containing 0.1% sodium azide.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.