SARS-ACSMS(C) (residues 1051-1076/1121-1154/1162-1190)

Ordering Information

Order Online »

JBS-PR-1105	100µg

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	-20°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	SARS-Associated Coronavirus Spike Mosaic S(C)
Application Notes	Recombinant SARS-ACSM Antigen may be used in ELISA and Western blots, excellent for detection of SARS with minimal specificity problems.
Formulation	liquid (Supplied in 25 mM Tris-HCl, 0.4% sarcosyl, 0,25% Triton X-100 and 50% glycerol)
MW	37 kDa
Purity	> 95 % (SDS-PAGE)
Source	recombinant, E. coli
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.