SEEBRIGHT® Green 496 dUTP (lyophilized)

Green 496 [5-Fluorescein] dUTP can replace TTP in reactions in which it serves as a substrate for *E. coli* DNA polymerase (holoenzyme and Klenow fragment), T4 and Taq DNA polymerases, reverse transcriptase (from AMV and M-MuLV) and terminal transferase. Fluorescently labeled probes can be prepared with this fluorescent nucleotide by a variety of methods including nick translation, random prime labeling, cDNA labeling and 3'-end labeling. Probes generated by these methods are suitable for use for the identification of specific sequences by *in situ* hybridization procedures on fixed cells and tissues by direct fluorescence detection. Green 496 dUTP can also be used for multicolor fluorescence labeling.

This labeled dUTP can be used with the Nick Translation DNA Labeling System 2.0 (Prod. No. ENZ-GEN111).

Citations: 8 View Online »

Ordering Information

Order Online »

ENZ-42831L-0050	50nmol

Manuals, SDS & CofA

View Online »

Handling & Storage

Use/Stability	As indicated on product label or CoA when stored as recommended. Stable for at least one year after receipt when stored as recommended.
Handling	Protect from light. Avoid freeze/thaw cycles.
Long Term Storage	-20°C
Shipping	Dry Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	5-Fluorescein dUTP
Appearance	Reddish orange solid.
Correction Factor (260nm)	0.32
Correction Factor (280nm)	0.2
Emission Maximum	520 nm
Excitation Maximum	496 nm
Extinction Coefficient	85,000 M-1 cm-1 (496 nm in TE [10 mM TRIS, pH 8.0, 1 mM EDTA])
Formulation	Lyophilized.
Purity	≥93% (HPLC)
Purity Detail	Purified by ion-exchange chromatography.
Quantity	Sufficient for approximately 98 reactions, following the recommended protocol of Prod. No. ENZ-GEN111.
Technical Info / Product Notes	Several of Enzo’s products and product applications are covered by US and foreign patents and patents pending.

Last modified: May 29, 2024

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.