Sp-2′-AHC-cAMPS-Agarose

PDE-resistant Sp-cAMPS immobilized on agarose gel

The PDE-resistant protein kinase A activator Sp-cAMPS immobilized on agarose by an aminohexylcarbamoyl spacer attached to the ribose 2' position of the ligand. The gel can be used for affinity chromatography of various cyclic nucleotide-responsive proteins, especially those which tolerate modification of the 2' ribose hydroxyl group such as some phosphodiesterases. This type of gel is also available with immobilized PDE-stable Rp-cAMPS (Rp-2'-AHC-cAMPS-Agarose (Prod. No. BLG-A049)) or with normal cAMP (2'-AHC-cAMP-Agarose (Prod. No. BLG-A058)).


Ordering Information

Order Online »

BLG-A050-06	0.6ml
BLG-A050-25	2.5ml
BLG-A050-60	6ml

Manuals, SDS & CofA

View Online »


NH2
N
N
N
N
O
O
S
P
O
O
O
Na+
H
N
O
N
H

Handling & Storage

Long Term Storage +4°C

Shipping Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name 2′-(6-Aminohexylcarbamoyl)adenosine-3′,5′-cyclic monophosphorothioate, Sp-isomer; immobilized on agarose gel

Couple Type Modified nucleotides

Technical Info / Product Notes For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.