Topoisomerase I Assay Kit

Assay kit contains reagents for routine detection of type I topoisomerase activity. DNA markers (included) allow unequivocal detection of the enzyme. Assay is specific for topoisomerase I; (eukaryotic type II topoisomerases require $MgCl_2$ and ATP). The standard kit is supplied with EDTA containing reaction buffer (minus Mg++) in order to prevent nucleolytic nicking of the supercoiled substrate (although we note that some eukaryotic type I activities are stimulated by Mg++). Assay will work even with crude extracts.

Ordering Information

Order Online »

TOP-TG1015-1	100 assays

Manuals, SDS & CofA

View Online »

Handling & Storage

Short Term Storage	+4°C
Long Term Storage	-20°C
Shipping	Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.