Topoisomerase IIα (human) Marker Set

The Western blot marker set for topoisomerase IIα is used to unambiguously locate the position of the 170 kDa topo II. The human 170 kDa is present in SDS-PAGE loading buffer.

Ordering Information

Order Online »

TOP-TG2011-3	20U

Manuals, SDS & CofA

View Online »

Handling & Storage

Handling	After opening, prepare aliquots and store at -80°C. Avoid freeze/thaw cycles.
Short Term Storage	-80°C
Long Term Storage	-80°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Contents	Topoisomerase IIα marker (20 units) in 1X SDS-PAGE load buffer. Marker contains the 170 kDa topoisomerase IIα.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.