Topoisomerase II Assay Kit

Assay kit is deigned to allow quick and specific assays for eukaryotic type II DNA topoisomerases. Facilitates the purification and characterization of type II enzymes and contains all reagents necessary for routine assays of type II topoisomerases.

Ordering Information

Order Online »

TOP-TG1001-2	250 assays

Manuals, SDS & CofA

View Online »

Handling & Storage

Short Term Storage +4°C

Long Term Storage -20°C

Shipping Ambient Temperature

Regulatory Status RUO - Research Use Only

Product Details

Contents

Kit Contents: (representative for 100 assay kit size)
-Kinetoplast DNA (KDNA) 20ug
-Marker linear KDNA in gel loading buffer
-Marker decatenated KDNA (topo II treated)
-10X Topoisomerase II assay buffer*
*NEW: To improve stability, the 10X Topo II Assay buffer is provided as two components: Buffer A
(no ATP) and Buffer B (ATP). The complete 10X Assay buffer is made fresh each time by combining
Buffer A and B.
-5X Stop buffer/ gel loading dye
-Detailed instruction manual

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.