Topoisomerase II Immunoprecipitation Kit

The immunoprecpitation (IP) kit contains the key reagents and protocol for immunoprecipitating p170 human topo IIα (Prod. No. TOP-2000H). The method should work well with other cells besides human cells. For example, mouse, rat and monkey cells will cross react with topo II antibody (Prod. No. TOP-2011-1). The antibody is specific for p170 (topo IIα) in human (HeLa) cells and although it is presumed that it will behave similarly in other species with cross reactive topo II isozymes, this should be carefully evaluated.

Ordering Information

Order Online »

TOP-TG1035	1Kit

Manuals, SDS & CofA

View Online »

Handling & Storage

Short Term Storage	-20°C
Long Term Storage	-20°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Quantity	>50 analyses

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.