Yeast DNA Preparation Kit (Spin column based), L pack

Spin column based genomic DNA purification from yeast

Ordering Information

Order Online »

JBS-PP-215L	5x50preparations

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage Ambient

Shipping Ambient Temperature

Regulatory Status

RUO - Research Use Only

Product Details

Application Notes

The obtained DNA is suitable for a variety of applications, including real-time PCR, southern blot analysis, genotyping and discovery or validation of SNP/SSR markers. Before start, add the following components (not included in the kit) as indicated on the respective bottle/tube: Lyticase Buffer to Lyticasedouble distilled water to RNase A and Proteinase K96-99 % Ethanol to the Washing Buffer

Technical Info / Product Notes

For the Original Manufacturer's data sheet please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.