ZytoDot® 2C SPEC USP6 Break Apart Probe

For the qualitative detection of translocations involving the human USP6 gene at 17p13.2by chromogenic in situhybridization (CISH)

The ZytoDot ® 2C SPEC USP6 Break Apart Probe is designed to detect translocations involving the chromosomal region 17p13.2 harboring the USP6 gene. Translocations affecting USP6 have been initially found in primary aneurysmal bone cysts, a benign, but locally aggressive bone lesion. USP6 rearrangements are restricted to spindle cells in primary ABC, indistinguishable from surrounding normal spindle cells.

Ordering Information

Order Online »

ZTV-C-3077-100	100µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	CISH
Quantity	Sufficient for 10 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.