ZytoDot® SPEC MDM2 Probe (Digoxigenin)

For the detection of the human MDM2 gene by chromogenic *in situ hybridization (CISH).*

The ZytoDot ® SPEC MDM2 Probe is designed for the detection of MDM2 gene amplifications found in more than 10% of human tumors. Due to the amplification of the respective chromosomal region, MDM2 is overexpressed in many human tumors such as soft tissue sarcomas, osteosarcomas, gliomas, NSCLC, gastric and breast carcinomas. Well-differentiated liposarcomas, the most common soft tissue tumors in adults, are characterized by the amplification of 12q-derived chromosomal material.

Ordering Information

Order Online »

ZTV-C-3012-400	400µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	Mouse Double Minute 2
Application	CISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 40 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.