ZytoFast® HPV Screening Probe (Digoxigenin)

For the detection of human Human Papilloma Virus (HPV) type

6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82

DNA by chromogenic *in situ* hybridization (CISH)

The ZytoFast HPV-CISH System is designed for the detection and discrimination of human papilloma virus (HPV) DNA in paraffin-embedded tissue sections or cell samples. Until now, approximately 10–15 HPV types are associated with lesions that can progress to cancer. The infection with the HPV hr-types can lead to development of cancer of the cervix, vulva, vagina, anus, or penis.

Ordering Information

Order Online »

ZTV-T-1144-400	400µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Handling	Protect from light.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	Human Papilloma Virus type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82
Application	CISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 40 tests, using 10 µl probe solution per test.
Specificity	HPV type 6/11/16/18/31/33/35/39/45/51/52/56/58/59/66/68/82
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.