# ZytoLight® FISH-Cytology Implementation Kit

**For fluorescence *in situ* hybridization (FISH) on cytology specimens using any Zyto*Light* FISH probe**

## Ordering Information

Order Online »

| | |
|---|---|
| **ZTV-Z-2099-20** | 20 tests |

## Manuals, SDS & CofA

View Online »

## Handling & Storage

| | |
|---|---|
| **Handling** | Protect from light. |
| **Long Term Storage** | +4°C |
| **Shipping** | Blue Ice |

## Regulatory Status

RUO - Research Use Only

## Product Details

| | |
|---|---|
| **Contents** | Cytology Pepsin Solution, 4 ml; 20x Wash Buffer TBS, 50 ml; 10x MgCl2, 50 ml; 10x PBS, 50 ml; Cytology Stringency Wash Buffer SSC, 500 ml; Cytology Wash Buffer SSC, 500 ml; DAPI/Antifade-Solution, 0.8 ml |
| **Technical Info / Product Notes** | For the Manual please click here. |

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.