ZytoLight® SPEC Angelman Dual Color Probe

For the qualitative detection of human UBE3A deletions and the human PML gene by fluorescence in situ hybridization (FISH)

The ZytoLight SPEC Angelmann Dual Color Probe (PL273) is intended to be used for the qualitative detection of human UBE3A gene deletions as well as the detection of the human PML gene in cytologic specimens by fluorescence in situ hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Cytology Implementation Kit (Prod. No. Z-2099-20).

Ordering Information

Order Online »

ZTV-Z-2319-50	50µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Handling	Protect from light.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	FISH
Quantity	Sufficient for 5 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.