ZytoLight® SPEC ERG Dual Color Break Apart Probe

For the detection of translocations involving the ERG gene at 21q22.2 by fluorescence *in situ* hybridization (FISH).

The ZytoLight® SPEC ERG Dual Color Break Apart Probe is designed to detect aberrations involving the ERG gene at 21q22.2 frequently detected in prostate cancers. ERG rearrangements have been observed in 40-60% of prostate cancers identified via prostate-specific antigen screening.

Ordering Information

Order Online »

ZTV-Z-2138-200	200µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	v-ETS Erythroblastosis Virus E26 Oncogene Homolog
Application	FISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 20 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.