ZytoLight® SPEC FGFR1 Dual Color Break Apart Probe

For the detection of translocations involving the FGFR1 gene at 8p11.22-p11.23 by fluorescence *in situ* hybridization (FISH).

The ZytoLight® SPEC FGFR1 Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal region 8p11.23-p11.22 harboring the FGFR1 gene. Translocations affecting FGFR1 are hallmarks of the 8p11 myeloproliferative syndrome also known as stem cell leukemia/lymphoma syndrome, an aggressive stem cell myeloproliferative neoplasm that is associated with eosinophilia, poor prognosis, T-cell lymphoma, and frequent progression to acute myeloid leukemia.

Ordering Information

Order Online »

ZTV-Z-2168-50	50µl
ZTV-Z-2168-200	200µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	Fibroblast Growth Factor Receptor 1, FGF Receptor 1
Application	FISH
Quantity	Sufficient for 5 tests (50 µl size) or 20 tests (200 µl size), using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.