ZytoLight® SPEC FOXO1 Dual Color Break Apart Probe

For the detection of translocations involving the FOXO1 gene at 13q14.11 by fluorescence *in situ* hybridization (FISH).

The ZytoLight® SPEC FOXO1 Dual Color Break Apart Probe is designed for the detection of specific translocations involving the chromosomal region 13q14.11 harboring the FOXO1 gene characteristic for alveolar rhabdomyosarcoma. Among solid tumors of the childhood, rhabdomyosarcoma is the most common soft tissue sarcoma. RMS are classified in two main categories: embryonal rhabdomyosarcoma and alveolar rhabdomyosarcoma. The alveolar histology is associated with a poorer prognosis.

Ordering Information

Order Online »

ZTV-Z-2139-50	50µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	Forkhead Box O1, FKHR
Application	FISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 5 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.