ZytoLight® SPEC PAX5 Dual Color Break Apart Probe

For the qualitative detection of translocations involving the human PAX5 gene at 9p13.2 by fluorescence *in situ* hybridization (FISH)

The ZytoLight® SPEC PAX5 Dual Color Break Apart Probe is designed to detect rearrangements involving the chromosomal region 9p13.2 harboring the PAX5 gene. PAX5 rearranged ALL patients were shown to respond well to treatment with prednisone. Hence, the identification of PAX5 rearrangements by FISH may be of therapeutic significance in ALL.

Ordering Information

Order Online »

ZTV-Z-2300-50	50µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	FISH
Quantity	Sufficient for 5 tests using 10 µl probe solution per test.
Technical Info / Product Notes	For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.