ZytoLight® SPEC SS18 Dual Color Break Apart Probe

For the detection of translocations involving SS18 gene at 18q11.2 by fluorescence *in situ* hybridization (FISH).

The ZytoLight® SPEC SS18 Dual Color Break Apart Probe is designed to detect translocations involving the chromosomal region 18q11.2 harboring the SS18 gene. Translocations involving the region 18q11.2 are found in over 90% of synovial sarcoma. Among soft tissue sarcomas, synovial sarcoma is one of the most common and classically occurs in the extremities of young adults even though, the occurrence of synovial sarcoma has also been described in a wide variety of anatomical locations and in all ages.

Ordering Information

Order Online »

ZTV-Z-2097-50	50µl
ZTV-Z-2097-200	200µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Handling	Protect from light.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	SYT
Application	FISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 5 tests (50 µl size) or 20 tests (200 µl size), using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.