ZytoLight® SPEC TFE3 Dual Color Break Apart Probe

For the detection of translocations involving the TFE3 gene at Xp11.23 by fluorescence *in situ* hybridization (FISH).

The ZytoLight® SPEC TFE3 Dual Color Break Apart Probe is designed to detect translocations involving the chromosomal region Xp11.23 harboring the TFE3 gene. Translocations involving the chromosomal region Xp11.2 are frequently detected in renal cell carcinomas (RCCs) which usually affect children and adolescents.

Ordering Information

Order Online »

ZTV-Z-2109-50	50µl
ZTV-Z-2109-200	200µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Alternative Name	Transcription Factor Binding To IGHM Enhancer 3, TFEA
Application	FISH
Formulation	Solution of salt (containing synthetic polynucleotides) in water
Quantity	Sufficient for 5 tests (50 µl size) or 20 tests (200 µl size), using 10 µl probe solution per test.
Technical Info / Product Notes	For the Original Manufacturer's data sheet please click here. For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.