ZytoLight® SPEC TP63/TBL1XR1 TriCheck Probe

For the qualitative detection of TBL1XR1-TP63 rearrangements by fluorescence in situ hybridization (FISH)

The ZytoLight SPEC TP63/TBL1XR1 TriCheck Probe (PL274) is intended to be used for the qualitative detection of rearrangements involving the human TP63 gene at 3q28 and the human TBL1XR1 gene at 3q26.32 in formalin-fixed, paraffin-embedded specimens by fluorescence in situ hybridization (FISH). The probe is intended to be used in combination with the ZytoLight FISH-Tissue Implementation Kit (Prod. No. Z-2028-5/-20).

Ordering Information

Order Online »

ZTV-Z-2320-50	50µl

Manuals, SDS & CofA

View Online »

Handling & Storage

Handling	Protect from light.
Long Term Storage	+4°C
Shipping	Blue Ice

Regulatory Status

RUO - Research Use Only

Product Details

Application	FISH
Quantity	Sufficient for 5 tests, using 10 µl probe solution per test.
Technical Info / Product Notes	For the Manual please click here.

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.