# ZytoLight® SPEC ZNF384 Dual Color Break Apart Probe

## For the qualitative detection of translocations involving the chromosomal region 12p13.31 harbouring the ZNF384 gene by fluorescence in situ hybridization (FISH)

The ZytoLight® SPEC ZNF384 Dual Color Break Apart Probe is designed to detect translocations involving the chromosomal region 12p13.31 harboring the ZNF384 gene. Rearrangements of the ZNF384 gene are recurrent in acute leukemia and are most frequently found in precursor B-cell acute lymphoblastic leukemia.

## Ordering Information

Order Online »

| | |
|---|---|
| **ZTV-Z-2275-50** | 50µl |

## Manuals, SDS & CofA

View Online »

## Handling & Storage

| | |
|---|---|
| **Long Term Storage** | +4°C |
| **Shipping** | Blue Ice |

## Regulatory Status

RUO - Research Use Only

## Product Details

| | |
|---|---|
| **Application** | FISH |
| **Quantity** | Sufficient for 5 tests, using 10 µl probe solution per test. |
| **Technical Info / Product Notes** | For the Manual please click here. |

ENZO LIFE SCIENCES, INC.
Phone: 800.942.0430
info-usa@enzolifesciences.com

European Sales Office
ENZO LIFE SCIENCES (ELS) AG
Phone: +41 61 926 8989
info-eu@enzolifesciences.com

Belgium, The Netherlands & Luxembourg
Phone: +32 3 466 0420
info-be@enzolifesciences.com

France
Phone: +33 472 440 655
info-fr@enzolifesciences.com

Germany
Phone: +49 7621 5500 526
info-de@enzolifesciences.com

UK & Ireland
Phone (UK customers): 0845 601 1488
Phone: +44 1392 825900
info-uk@enzolifesciences.com

For local distributors and detailed product information visit us online: enzo.com All product names, logos, and brands are the property of their respective owners.